

Methods for Lab Data: C-Reactive Protein: Exam 6 Offspring

Analyte measured	CRP
Cohort	offspring
Exam Number	6
Sample Type	frozen serum
Storage Temp	-80C
Method Type	Particle Enhanced Immunonephelometry
Instrumentation (manufacturer/model)	Dade Behring BN 100
Kit/Reagent Name	Dade Behring - high sensitivity CRP reagent
Kit/Reagent Manufacturer	Dade
Measuring Range	0.16 - 550 mg/L
Run in single/duplicate/other	single
Inter-assay CV	3.4%
Intra-assay CV	3.8%
Run dates	2/04 - 3/04