

Methods for Lab Data: Urine Creatinine, Albumin and Sodium: Exam 6 Offspring

Analyte measured	urine creatinine
Cohort	offspring
Exam Number	6
Sample Type	random urine
Storage Temp	-20C
Method Type	colorimetric
Instrumentation (manufacturer/model)	Roche Hitachi 911
Kit/Reagent Name	Creat Rate-Blanked
Kit/Reagent Manufacturer	BM Roche
Measuring Range	
Run in single/duplicate/other	unknown
Inter-assay CV	
Intra-assay CV	
Run dates	12/98 - 1/99

Analyte measured	urine albumin
Cohort	offspring
Exam Number	6
Sample Type	random urine
Storage Temp	-20C
Method Type	turbidometric
Instrumentation (manufacturer/model)	Roche-Hitachi 911
Kit/Reagent Name	Albumin
Kit/Reagent Manufacturer	BM Roche
Measuring Range	
Run in single/duplicate/other	unknown
Inter-assay CV	
Intra-assay CV	
Run dates	12/98 - 1/99

Analyte measured	urine sodium
Cohort	offspring
Exam Number	6
Sample Type	random urine
Storage Temp	-20C
Method Type	ISE
Instrumentation (manufacturer/model)	Roche-Hitachi 911
Kit/Reagent Name	ISE NA,K,CL
Kit/Reagent Manufacturer	BM Roche
Measuring Range	
Run in single/duplicate/other	unknown
Inter-assay CV	
Intra-assay CV	
Run dates	12/98 - 1/99